B-CURE® LASER **SPORT** INSTRUCTIONS FOR USE

GOOD
ENERGIES®

Content


English	Instructions for use.....	03
Français	Mode d'emploi	11
Deutsch	Gebrauchsanweisung	19
Italiano	Istruzioni per l'uso.....	27
Русский	Инструкция по эксплуатации	35
Español	Instrucciones de uso	43
Português	Instruções de utilização.....	51

For better comfort it is recommended to use
B-Cure adjustable stand (not included)

Content:

1. Component legend
2. Battery charging
3. Battery warnings
4. Safety precautions
5. Contra indications
6. Operating instructions
7. Treatment protocol
8. Control panel and troubleshooting
9. Technical specifications
10. Device maintenance
11. Warranty

1. Component legend**Package content:**

- B-Cure Laser device
- B-Cure case
- Instructions For Use
- Batteries
- Charger
- Adaptors

Upon opening the package, please verify that it contains all the above mentioned equipment. Please make sure the device is unharmed. If anything is missing or damaged, please contact the distributor immediately after purchasing.

2. Battery charging

Inside the battery compartment there are 3 (type AAA) batteries partly charged. Complete the charging by plugging the charger into an electrical outlet and connect the charger to the socket. Battery charge time is approximately 7 hours, or until the word "FULL" appears on the screen. We recommend you charge the device overnight once every three days, or when the display blinks/the device stops functioning.

**3. Battery warnings**

- The AC/DC charger that comes with the device is the only charger allowed for use.
- Do not take apart the batteries, heat them over 100°C or burn them.
- Discard the batteries at an appropriate recycling center.
- Do not turn the device on while it is being charged.
- Do not, under any circumstances, charge ordinary, non-rechargeable batteries; it might cause a fire.

4. Safety precautions

- ⚠ • Never aim the laser beam at the eyes.
• It is always recommended to consult your doctor.
• If you are pregnant do not use the device to treat the stomach.
• Do not treat an area with a malignant tumor.

5. Contraindications

None

6. Operating instructions

6.1 Switch the device ON

Press the main power switch (No. 1). After a short beep the display will switch on and a timer will be displayed (No. 3).

6.2 Set treatment time

Press the treatment time setting button (No. 2): on its right side to increase time (+) or its left side to decrease time (-). The timer is adjustable in increments of 30 seconds.

6.3 Activate the laser beam and begin treatment

Step 1: press the “treatment start push-button” (no. 4 in the component legend illustration) with one hand, and make sure the green light at the front of the device is on. Step 2: continue pressing this button and simultaneously

press, using an additional finger on the pressure switch area (no. 1 in the illustration). Step 3: release both buttons at the same time. Now the laser beam is active and treatment can begin. Place the device gently in full contact with the skin. Do not treat on top of clothes or bandaging. To stop the functioning of the device, briefly press the “treatment start push-button” (no. 4 in the illustration).

During treatment, you will see a green light indicating the device is functioning. “Beep” sounds at 3 second intervals indicate that the device is functioning and that the stopwatch is counting down until reaching “00:00” at which time the device’s functioning will automatically end and the “Beep” sound will stop. Hold the device on the skin steadily and motionless during the whole duration of the treatment.

Demonstration of the activation of the laser beam after setting treatment time

After setting the treatment time, press the “treatment start push-button”

Continue pressing and simultaneously press using an additional finger on the pressure switch area

Release both buttons simultaneously. Now the green light is on and treatment can begin

6.4 End of treatment, reset and turn OFF

To switch off the device, press the main power switch until a beep sounds. Release the button and the display will go blank. If the planned treatment time has not been fully utilized, the timer will display the remaining time at startup. To reset the timer, switch off the device, then press the main power switch for longer than 2 seconds while switching on the device.

7. Treatment protocol

7.1 Intended use

B-Cure Laser is intended for the management of pain including: Chronic back pain; Chronic neck pain; Knee pain; Carpal Tunnel Syndrome (CTS); Fibromyalgia; Myofascial Pain Syndrome (MPS); Temporomandibular Joint Disorder; May help as an adjunctive treatment for acute and chronic wounds.

7.2 Treatment course and dosage

Start treatment in 1.5 minute dosage for each treatment point and increase the dosage gradually within 24 hours until you reach 8 min. Treat twice a day, consistently. If you experience increased pain following treatment, stop the treatment for 24 hours then start again with a reduced dosage. You may increase the dosage to 4 treatments a day. In cases of acute pain, trauma or non-chronic aches, you may apply a dosage of 8 minutes from the first treatment.

7.3 Coverage

The laser beam size is 4.5 cm long and 1 cm wide. If the area that requires treatment is larger than the 4.5/1 cm, repeat treatments on the adjacent area until you've covered the whole area.


7.4 Treatment recommendations for common conditions

The recommended time for each point of treatment is 8 minutes, 2-4 times a day.

For all recommended treatment therapies, consult the human body illustration, page 7.

- 1. Chronic back pain** - Treat lengthwise along the painful vertebrae of the lower back, to the left and the right of the spine. For Ischia - pain radiating to the lower limb - Identify the focal point of the pain by touch and treat it.
- 2. Chronic neck pain** - Treat **a)** Lengthwise, along the vertebra. **b)** Treat each side of the neck, lengthwise; under the ear lobe, and exactly at the indentation behind the jaw.
- 3. Knee pain** - Treat the sore knee **a)** On the back of the knee along the inner crease. **b)** On the outer side of the knee, lengthwise. **c)** On the inner side of the knee, lengthwise. You can also treat the front of the knee, under the kneecap.
- 4. Carpal tunnel syndrome (CTS)** - Treat the sore wrist **a)** On the inside, along the width of the joint - on the crease. **b)** On the inside, perpendicular to 'a', lengthwise. **c)** Treat the length of the thumb, along the pad. You can also treat the back of the hand, at the gap between the thumb and the bone of the index finger towards the joint, and from both sides of the joint, lengthwise.
- 5. Fibromyalgia** - Treat the aching area.
- 6. Myofascial pain syndrome (MPS)** - Treat the aching area.
- 7. Temporomandibular Joint Disorder** - Treat along the length and width of the TemporoMandibular joint.
- 8. Adjunctive treatment for acute and chronic wounds** - Treat over the entire wound area.


Tip: Always start from the most painful point and then move to other spots, so as to ensure you are treating the central problem first, to the best of your knowledge.

The size of the laser beam of the
B-CURE® LASER

Human body illustration with recommended treatment points:

- Note the illustrated beam direction
- You can treat any point in the body except the eyes
- In case of a raw wound or skin tissue, do not place the device in full contact with the skin. In any other case, place the device on the skin itself.
- Do not treat on top of clothes or bandaging.
- Do not move the device during treatment.
- It is advisable to treat all the recommended points, to ensure treatment efficiency.

8. Control panel and troubleshooting

Display Status

Batteries charging

Batteries fully charged

The icon-* is displayed to indicate the device is currently operating

Full battery power

Low battery power, two beeps will sound after startup as a reminder to recharge. Charge the batteries

No battery power remains. Charge the batteries immediately

The device will shut down automatically after 10 beeps due to lack of battery power. Charge the batteries

Batteries are not charging. Switch the device off then on again to fix this error.

If the problem persists, switch off the device, remove the batteries and then reinsert them.

If the problem persists, ship or deliver personally the device to the service center.

The indication "LOW" is displayed when laser power level is below the required level.

Treatment should be stopped. Cause may be one of the following:

1. Battery power is too low.
2. Diode lifespan has reached its end and should be changed at the service center.
3. The device is operated in a different environment than for which it was designed, e.g. temperatures higher than 40°C. Charge the batteries, switch on the device and see if the writing "LOW" has disappeared. If not, verify that room temperature is not exceptionally high, and if necessary move away from the hot area. If not higher than 40°C and the writing "LOW" is still displayed - ship the device to the service center for diode change.

Additional signals for errors of this type:

Hardware/software error has occurred. Switch the device off, then on again, to fix this error.

If the problem persists, switch off the device, remove the batteries and then reinsert them.

If the problem persists, ship or deliver personally the device to the service center.

9. Technical specifications

Property	Units	Description, value	Notes
Laser type		GaAlAs	Solid state laser diode
Maximum power	mW	250	Max
Wave length (infrared)	nM	808	
Pulse frequency	kHz	14	
Laser pulse duration	μs	22	
Diode lifespan	hours	Over 3000	
Energy per minute	Joul	4.5	
Laser beam dimensions (LxH)	mm	Length 45 x width 10	at opening
Beam divergence	degrees	±3	
Green control light indicator	nM	520	Green light beam signaling device is operating
Timer setting	minutes	00:00 to 99:30	
Timer intervals	seconds	30	Increment or decrement
LCD display dimensions	mm	26 x 15	Visible area
Timer display		4 digits (88:88)	1-second intervals
Beep warning		yes	Display: Device on, off, operating, low laser power, errors
Power supply	W	3.2	maximum
Automatic shutdown		2 options	automatic shutdown when not in use for 10 minutes; 5 minutes when timer is set to 0 (00:00)
Battery			3 AAA rechargeable batteries (NiMH)
Batteries lifespan	hours	over 2	with fully charged batteries
Battery charging time	hours	5-7	with discharged batteries
Batteries charging current	mA	300	minimum
Power supply	VAC	100 to 240	AC/DC charging adapter 115/230 VAC, 65/50 Hz Supply 9V, 500mA

9. Technical specifications - continued

Property	Units	Description, value	Notes
Dimensions	mm	200 x 70 x 40	
Charger output voltage	VDC	9	
Weight	grams	175	
Operating conditions	C°, %	temp. 0° to 40°, humidity 15% to 90%	without condensation
Storage conditions	C°, %	temp. -20° to 55°, humidity 10% to 95%	without condensation

10. Device maintenance

1. Clean the laser beam aperture (No. 5), using a 70% alcohol solution . Clean between sessions if needed and before handing it to another person's usage.
2. While cleaning, be sure to prevent wetting the lens.
3. Keep the device in proper storage and work conditions as specified in the technical specification table, do not leave the device in the car in hot weather.

11. Warranty

The warranty is for 12 months according to company policy, as is detailed on warranty certificate. If batteries and/or a charger are supplied with the device, they are not included in the warranty. The manufacturer/ distributor is responsible for performance, safety and reliability of the equipment only on condition that all repairs are carried out exclusively by Company - authorized personnel.

Consult accompanying documents 	Temperature limitation  0°C - 40°C	Humidity limitation  10% - 95%
Electrical and electronic equipment 	CE marking of conformity 	Caution, consult accompanying documents 
Type BF 	IP22 – International Protection	Manufacturer 
Recycling 	Authorized representative in the European community	EC REP
3 x 1.2V AAA Rechargeable batteries 9V, 300 mA (For charging via AC/ DC adaptor)		

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Manufacturer: Biocare Enterprise LTD
Distributor: Erika Carmel LTD